

CHRND Mouse mAb

Catalog No: #63477

Package Size: #63477-1 50ul #63477-2 100ul


Description

Product Name	CHRND Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CHRND (AA: extra 22-245) expressed in E. Coli.
Target Name	CHRND
Other Names	ACHRD; CMS2A; CMS3A; CMS3B; CMS3C; FCCMS; SCCMS
Accession No.	Q07001
Calculated MW	58.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CHRND mouse mAb against C6 (1) cell lysate.



Flow cytometric analysis of SK-N-SH cells using CHRND mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.