

CHRNE Mouse mAb

Catalog No: #63478

Package Size: #63478-1 50ul #63478-2 100ul


Description

Product Name	CHRNE Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human CHRNE (AA: extra 21-239) expressed in E. Coli.
Target Name	CHRNE
Other Names	ACHRE; CMS1D; CMS1E; CMS2A; CMS4A; CMS4B; CMS4C; FCCMS; SCCMS
Accession No.	Q04844
Calculated MW	54.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CHRNE mouse mAb against C6 (1) cell lysate.



Flow cytometric analysis of SK-N-SH cells using CHRNE mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.