

CIB1 Mouse mAb

Catalog No: #63480

Package Size: #63480-1 50ul #63480-2 100ul


Description

Product Name	CIB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of CIB1 expressed in E. Coli.
Target Name	CIB1
Other Names	CIB; KIP; KIP1; SIP2-28
Accession No.	Q99828
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human thalamus (left) and glioma (right) tissue, showing membrane localization using CIB1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.