

CIDEC Mouse mAb

Catalog No: #63481

Package Size: #63481-1 50ul #63481-2 100ul


Description

Product Name	CIDEC Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CIDEC (AA: 53-141) expressed in E. Coli.
Target Name	CIDEC
Other Names	CIDE3; FPLD5; FSP27; CIDE-3
Accession No.	Q96AQ7
Calculated MW	26.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using CIDEC mouse mAb against HEK293 (1), A431 (2), and HCT116 (3) cell lysate.



Immunofluorescence analysis of HeLa cells using CIDEc mouse mAb (green). Blue: DAPI fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using CIDEc mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.