

CLDN6 Mouse mAb

Catalog No: #63487

Package Size: #63487-1 50ul #63487-2 100ul


Description

Product Name	CLDN6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CLDN6 (AA: extra(29-81) and extra(138-160)) expressed in E. Coli.
Target Name	CLDN6
Other Names	CLDN6
Accession No.	P56747
Calculated MW	23.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200-1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Hepg2 cells using CLDN6 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.