

CLEC5A Mouse mAb

Catalog No: #63489

Package Size: #63489-1 50ul #63489-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLEC5A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CLEC5A (AA: extra(28-188)) expressed in E. Coli.
Target Name	CLEC5A
Other Names	MDL1; MDL-1; CLECSF5
Accession No.	Q9NY25
Calculated MW	21.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CLEC5A mouse mAb against THP-1 (1) cell lysate.



Flow cytometric analysis of U937 cells using CLEC5A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.