

CNN1 Mouse mAb

Catalog No: #63492

Package Size: #63492-1 50ul #63492-2 100ul


Description

Product Name	CNN1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse,Rat,Monkey,Rabbit
Immunogen Description	Purified recombinant fragment of human CNN1 (AA: 16-165) expressed in E. Coli.
Target Name	CNN1
Other Names	SMCC; Sm-Calp
Accession No.	P51911
Calculated MW	33.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/100 - 1/500ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using CNN1 mouse mAb against C2C12(1) COS-7 (2) HELA(3) cell lysate.




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CNN1 mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using CNN1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Mouse spleen(A) Mouse muscle(B) using CNN1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat spleen(A) Rat muscle(B) using CNN1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit lung using CNN1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.