

COL2A1 Mouse mAb

Catalog No: #63495

Package Size: #63495-1 50ul #63495-2 100ul


Description

Product Name	COL2A1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human COL2A1 (AA: 1222-1487) expressed in mammalian.
Target Name	COL2A1
Other Names	AOM; ANFH; SEDC; STL1; COL11A3
Accession No.	P02458
Calculated MW	142kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using COL2A1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using COL2A1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using COL2A1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.