

COTL1 Mouse mAb

Catalog No: #63497

Package Size: #63497-1 50ul #63497-2 100ul


Description

Product Name	COTL1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human COTL1 expressed in E. Coli.
Target Name	COTL1
Other Names	CLP
Accession No.	Q14019
Calculated MW	16kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human breast cancer tissues using COTL1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using COTL1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.