

COX4I1 Mouse mAb

Catalog No: #63499

Package Size: #63499-1 50ul #63499-2 100ul


Description

Product Name	COX4I1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human COX4I1 expressed in E. Coli.
Target Name	COX4I1
Other Names	COX4; COXIV; COX4-1; MGC72016; COX4I1
Accession No.	P13073
Calculated MW	19kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Western blot analysis using COX4I1 mouse mAb against HEK293 (1), A549 (2) and PC12 (3) cell lysate.



Immunofluorescence analysis of PANC-1 cells using COX4I1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of K562 cells using COX411 mouse mAb (blue) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.