

CRKL Mouse mAb

Catalog No: #63504

Package Size: #63504-1 50ul #63504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CRKL Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CRKL expressed in E. Coli.
Target Name	CRKL
Other Names	CRKL
Accession No.	P46109
Calculated MW	39kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Flow cytometric analysis of NIH/3T3 cells using CRKL mouse mAb (blue) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.