

CTCF Mouse mAb

Catalog No: #63517

Package Size: #63517-1 50ul #63517-2 100ul


Description

Product Name	CTCF Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human CTCF expressed in E. Coli.
Target Name	CTCF
Other Names	CTCF
Accession No.	P49711
Calculated MW	83kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using CTCF mouse mAb against A31 (1), MCF-7 (2), Hela (3), HCT116 (4), Jurkat (5), NIH/3T3 (6), and Cos7 (7) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.