

CTTN Mouse mAb

Catalog No: #63521

Package Size: #63521-1 50ul #63521-2 100ul


Description

Product Name	CTTN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human CTTN expressed in E. Coli.
Target Name	CTTN
Other Names	EMS1; FLJ34459; CTTN
Accession No.	Q14247
Calculated MW	80kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using CTTN mouse mAb against Hela (1), A431 (2), MCF-7 (3), SR-BR-3 (4), HepG2 (5) and NIH/3T3 (6) cell lysate.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues (left) and tonsil tissues (right) using CTTN mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.