

DAPK3 Mouse mAb

Catalog No: #63529

Package Size: #63529-1 50ul #63529-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DAPK3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DAPK3 (AA: 28-161) expressed in E. Coli.
Target Name	DAPK3
Other Names	ZIP; ZIPK
Accession No.	O43293
Calculated MW	52.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using DAPK3 mouse mAb against A431 (1), K562 (2) cell lysate.



Flow cytometric analysis of A431 cells using DAPK3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using DAPK3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.