

DCLK2 Mouse mAb

Catalog No: #63532

Package Size: #63532-1 50ul #63532-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCLK2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DCLK2 (AA: 652-766) expressed in mammalian.
Target Name	DCLK2
Other Names	CL2; DCK2; CLIK2; DCDC3; CLICK2; DCDC3B; DCAMKL2; CLICK-II
Accession No.	Q8N568
Calculated MW	83.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of U937 cells using DCLK2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human brain tissues using DCLK2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.