

DCTN4 Mouse mAb

Catalog No: #63534

Package Size: #63534-1 50ul #63534-2 100ul


Description

Product Name	DCTN4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DCTN4 (AA: 57-298) expressed in E. Coli.
Target Name	DCTN4
Other Names	P62; DYN4
Accession No.	Q9UJW0
Calculated MW	52.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Western blot analysis using DCTN4 mouse mAb against Raw264.7 (1) and NIH3T3 (2) cell lysate.



Immunofluorescence analysis of HepG2 cells using DCTN4 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of A549 cells using DCTN4 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using DCTN4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.