

DDX20 Mouse mAb

Catalog No: #63538

Package Size: #63538-1 50ul #63538-2 100ul


Description

Product Name	DDX20 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human DDX20 (AA: 725-824) expressed in E. Coli.
Target Name	DDX20
Other Names	DP103; GEMIN3
Accession No.	Q9UHI6
Calculated MW	92.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using DDX20 mouse mAb against Jurkat (1), Hela (2), NTERA-2 (3), HL7702 (4), K562 (5), and C6 (6) cell lysate.



Immunofluorescence analysis of HeLa cells using DDX20 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using DDX20 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.