

DIS3L2 Mouse mAb

Catalog No: #63547

Package Size: #63547-1 50ul #63547-2 100ul


Description

Product Name	DIS3L2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DIS3L2 (AA: 27-250) expressed in E. Coli.
Target Name	DIS3L2
Other Names	FAM6A; PRLMNS
Accession No.	Q8IYB7
Calculated MW	99.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using DIS3L2 mouse mAb against HeLa (1) and HepG2 (2) cell lysate.



Flow cytometric analysis of Jurkat cells using DIS3L2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using DIS3L2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using DIS3L2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.