

DLG4 Mouse mAb

Catalog No: #63550

Package Size: #63550-1 50ul #63550-2 100ul


Description

Product Name	DLG4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4549
Isotype	Mouse IgG1
Applications	WB,IHC,IF,FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human DLG4 (AA: 54-300) expressed in HEK293-6e cells supernatant.
Target Name	DLG4
Other Names	MRD62; PSD95; SAP90; SAP-90
Accession No.	P78352
Calculated MW	80.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using DLG4 mouse mAb against Mouse brain (1), and Rat brain (2) tissue lysate.



Flow cytometric analysis of LNCAP cells using DLG4 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human brain tissues using DLG4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.