

DNM1L Mouse mAb

Catalog No: #63558

Package Size: #63558-1 50ul #63558-2 100ul


Description

Product Name	DNM1L Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human DNM1L (AA: 69-213) expressed in E. Coli.
Target Name	DNM1L
Other Names	DLP1; DRP1; DVLP; EMPF; VPS1; DYMPLE; HDYNIV; DYNIV-11
Accession No.	O00429
Calculated MW	81.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using DNM1L mouse mAb against PC-12 (1), LNCAP (2) and NIH/3T3 (3) cell lysate.



Immunofluorescence analysis of HeLa cells using DNM1L mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of HEK293 cells using DNM1L mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.