

DNMT3B Mouse mAb

Catalog No: #63561

Package Size: #63561-1 50ul #63561-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DNMT3B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DNMT3B (AA: 1-150) expressed in E. Coli.
Target Name	DNMT3B
Other Names	ICF; ICF1; M.HsallIB
Accession No.	Q9UBC3
Calculated MW	95.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of K562 cells using DNMT3B mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using DNMT3B mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using DNMT3B mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.