

DNMT3L Mouse mAb

Catalog No: #63562

Package Size: #63562-1 50ul #63562-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DNMT3L Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DNMT3L (AA: 147-386) expressed in E. Coli.
Target Name	DNMT3L
Accession No.	Q9UJW3
Calculated MW	43.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using DNMT3L mouse mAb (green) and negative control (red).



Flow cytometric analysis of HepG2 cells using DNMT3L mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.