

DSG3 Mouse mAb

Catalog No: #63565

Package Size: #63565-1 50ul #63565-2 100ul


Description

Product Name	DSG3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DSG3 (AA: 55-159) expressed in E. Coli.
Target Name	DSG3
Other Names	PVA; CDHF6
Accession No.	P32926
Calculated MW	107.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using DSG3 mouse mAb against A431 cell lysate.



Flow cytometric analysis of A431 cells using DSG3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using DSG3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal tissues using DSG3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.