

DYRK4 Mouse mAb

Catalog No: #63569

Package Size: #63569-1 50ul #63569-2 100ul


Description

Product Name	DYRK4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human DYRK4 (AA: 303-515) expressed in E. Coli.
Target Name	DYRK4
Accession No.	Q9NR20
Calculated MW	59.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using DYRK4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human pancreas tissues using DYRK4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.