

E7 Mouse mAb

Catalog No: #63572

Package Size: #63572-1 50ul #63572-2 100ul


Description

Product Name	E7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Virus
Immunogen Description	Purified recombinant fragment of E7 expressed in E. Coli.
Target Name	E7
Other Names	E7; HpV16gp2
Accession No.	P03129
Calculated MW	21kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunohistochemical analysis of paraffin-embedded human placenta tissues (left) and ovarian cancer (right) using E7 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.