

EGF Mouse mAb

Catalog No: #63575

Package Size: #63575-1 50ul #63575-2 100ul


Description

Product Name	EGF Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human EGF expressed in E. Coli.
Target Name	EGF
Other Names	URG; HOMG4; EGF
Accession No.	P01133
Calculated MW	134kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Immunohistochemical analysis of paraffin-embedded human breast cancer, Lung breast tissues using EGF mouse mAb



Immunohistochemical analysis of paraffin-embedded human Uterus tissues using anFtEG mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.