

EhpB1 Mouse mAb

Catalog No: #63578

Package Size: #63578-1 50ul #63578-2 100ul


Description

Product Name	EhpB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphB1 (aa19-133) expressed in E. Coli.
Target Name	EhpB1
Other Names	ELK; NET; Hek6; EPHT2; EPHB1
Accession No.	P54762
Calculated MW	110kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using EhpB1 mouse mAb against MDA-MB-468 (1), MDA-MB-453 (2), MCF-7 (3), T47D (4) and SKBR-3 (5) cell lysate.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.