

ELANE Mouse mAb

Catalog No: #63588

Package Size: #63588-1 50ul #63588-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ELANE Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human ELANE (AA: 140-267) expressed in E. Coli.
Target Name	ELANE
Other Names	GE; NE; HLE; HNE; ELA2; SCN1; PMN-E
Accession No.	P08246
Calculated MW	28.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ELANE mouse mAb against U937 (1), SPC-A-1 (2), and COS7 (3) cell lysate.



Flow cytometric analysis of Hela cells using ELANE mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using ELANE mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.