

EphA1 Mouse mAb

Catalog No: #63591

Package Size: #63591-1 50ul #63591-2 100ul


Description

Product Name	EphA1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphA1 expressed in E. Coli.
Target Name	EphA1
Other Names	EPH; EPHT; EPHT1
Accession No.	P21709
Calculated MW	108kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using EphA1 mouse mAb against A549 (1) and Hela (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovary carcinoma (left) and breast carcinoma (right), showing cytoplasmic localization using EphA1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.