

EphA2 Mouse mAb

Catalog No: #63592

Package Size: #63592-1 50ul #63592-2 100ul


Description

Product Name	EphA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgM
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphA2 expressed in E. Coli.
Target Name	EphA2
Other Names	ECK; EPHA2
Accession No.	P29317
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using EphA2 mouse mAb against NIH/3T3 cell lysate.



Immunohistochemical analysis of paraffin-embedded human skin carcinoma (left) and pancreas carcinoma (right) tissue, showing cytoplasmic localization using EphA2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.