

EphA7 Mouse mAb

Catalog No: #63593

Package Size: #63593-1 50ul #63593-2 100ul


Description

Product Name	EphA7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphA7 (aa27-210) expressed in E. Coli.
Target Name	EphA7
Other Names	EHK3; HEK11
Accession No.	Q15375
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human Breast tissues using EPHA7 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.