

EphB3 Mouse mAb

Catalog No: #63595

Package Size: #63595-1 50ul #63595-2 100ul


Description

Product Name	EphB3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphB3 (aa39-212) expressed in E. Coli.
Target Name	EphB3
Other Names	ETK2; HEK2; TYRO6
Accession No.	P54753
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung squamous cell carcinoma (A), lung adenocarcinoma (B), colon carcinoma (C), breast carcinoma (D), normal sublingual gland (E), normal rectal (F), showing membrane localization with DAB staining using EphB3 mouse mAb.



Immunohistochemical analysis of paraffin-embedded human Adrenal tissues using EPHB3 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.