

EphB4 Mouse mAb

Catalog No: #63597

Package Size: #63597-1 50ul #63597-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EphB4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of EphB4 (aa562-612) expressed in E. Coli.
Target Name	EphB4
Other Names	HTK; MYK1; TYRO11
Accession No.	P54760
Calculated MW	108kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using EphB4 mouse mAb against Jurkat (1) and HEK293 (2) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.