

EPHB4 Mouse mAb

Catalog No: #63598

Package Size: #63598-1 50ul #63598-2 100ul


Description

Product Name	EPHB4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant extracellular fragment of human EPHB4 fused with hlgGFc tag expressed in HEK293 cell line.
Target Name	EPHB4
Other Names	HTK; MYK1; TYRO11
Accession No.	P54760
Calculated MW	108kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ICC: 1/200 - 1/1000 ELISA: 1/10000

Images



Confocal Immunofluorescence analysis of methanol-fixed HEK293 cells trasfected with EphB4-hlgGFc using EphB4 mouse mAb (green), showing membrane localization. Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.