

EPO Mouse mAb

Catalog No: #63600

Package Size: #63600-1 50ul #63600-2 100ul


Description

Product Name	EPO Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human EPO expressed in E. Coli.
Target Name	EPO
Other Names	EP; MVCD2; MGC138142; EPO
Accession No.	P01588
Calculated MW	21kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Immunofluorescence analysis of GC7901 cells using EPO mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.