

ER Mouse mAb

Catalog No: #63601

Package Size: #63601-1 50ul #63601-2 100ul

Description

Product Name	ER Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of ER expressed in E. Coli.
Target Name	ER
Other Names	ER; ESR; Era; ESRA; NR3A1; DKFZp686N23123; ESR1
Accession No.	P03372
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.