

ERBB3 Mouse mAb

Catalog No: #63603

Package Size: #63603-1 50ul #63603-2 100ul


Description

Product Name	ERBB3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of ERBB3 (aa1175-1275) expressed in E. Coli.
Target Name	ERBB3
Other Names	HER3; LCCS2
Accession No.	P21860
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human cerebra (left) and breast carcinoma tissue (right), showing cytoplasmic and membrane location with DAB staining using ERBB3 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.