

ERK3 Mouse mAb

Catalog No: #63606

Package Size: #63606-1 50ul #63606-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ERK3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ERK3 expressed in E. Coli.
Target Name	ERK3
Other Names	ERK3; PRKM6; p97MAPK; HsT17250; DKFZp686F03189; MAPK6
Accession No.	Q16659
Calculated MW	105kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Flow cytometric analysis of Hela cells using ERK3 mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.