

ESRRA Mouse mAb

Catalog No: #63609

Package Size: #63609-1 50ul #63609-2 100ul


Description

Product Name	ESRRA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human ESRRA (AA: 198-376) expressed in E. Coli.
Target Name	ESRRA
Other Names	ERR1; ERRa; ESRL1; NR3B1; ERRalpha
Accession No.	P11474
Calculated MW	45.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ESRRA mouse mAb against C6 (1) cell lysate.



Flow cytometric analysis of Hela cells using ESRRA mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.