

ETV5 Mouse mAb

Catalog No: #63613

Package Size: #63613-1 50ul #63613-2 100ul


Description

Product Name	ETV5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ETV5 expressed in E. Coli.
Target Name	ETV5
Other Names	ERM; ETV5
Accession No.	P41161
Calculated MW	58kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using ETV5 mouse mAb against Jurkat (1), NIH/3T3 (2) and MCF-7 (3) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.