

EZH2 Mouse mAb

Catalog No: #63615

Package Size: #63615-1 50ul #63615-2 100ul


Description

Product Name	EZH2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human EZH2 (AA: 1-194) expressed in E. Coli.
Target Name	EZH2
Other Names	WVS; ENX1; EZH1; KMT6; WVS2; ENX-1; EZH2b; KMT6A
Accession No.	Q15910
Calculated MW	85.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using EZH2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using EZH2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using EZH2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.