

F2R Mouse mAb

Catalog No: #63617

Package Size: #63617-1 50ul #63617-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	F2R Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human F2R (AA: 42-176) expressed in E. Coli.
Target Name	F2R
Other Names	TR; HTR; CF2R; PAR1; PAR-1
Accession No.	P25116
Calculated MW	47.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using F2R mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.