

FAK Mouse mAb

Catalog No: #63620

Package Size: #63620-1 50ul #63620-2 100ul


Description

Product Name	FAK Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human FAK expressed in E. Coli.
Target Name	FAK
Other Names	FAK; FADK; FAK1; FRNK; pp125FAK; PTK2
Accession No.	Q05397
Calculated MW	119kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using FAK mouse mAb against A549 (1) and NIH/3T3 (2) cell lysate.



Immunofluorescence analysis of B16 cells using FAK mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human cervixes tumour using FAK mouse mAb with DAB staining

Note: This product is for in vitro research use only and is not intended for use in humans or animals.