

FAK Mouse mAb

Catalog No: #63621

Package Size: #63621-1 50ul #63621-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FAK Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FAK expressed in E. Coli.
Target Name	FAK
Other Names	FAK; FADK; FAK1; FRNK; pp125FAK; PTK2
Accession No.	Q05397
Calculated MW	119kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images



Immunofluorescence analysis of A549 cells using FAK mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of Raji cells using FAK mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.