

FCER1A Mouse mAb

Catalog No: #63624

Package Size: #63624-1 50ul #63624-2 100ul


Description

Product Name	FCER1A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FCER1A (AA: 42-103) expressed in E. Coli.
Target Name	FCER1A
Other Names	FCE1A; FcERI
Accession No.	P12319
Calculated MW	29.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using FCER1A mouse mAb against SW620 (1), A549 (2), and A431 (3) cell lysate.



Flow cytometric analysis of HEK293 cells using FCER1A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.