

FGF2 Mouse mAb

Catalog No: #63629

Package Size: #63629-1 50ul #63629-2 100ul


Description

Product Name	FGF2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of FGF2 expressed in E. Coli.
Target Name	FGF2
Other Names	BFGF; FGFB; HBGF-2; FGF2
Accession No.	P09038
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human rectum adenocarcinoma tissue showing cytoplasmic localization using FGF2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.