

FGF2 Mouse mAb

Catalog No: #63631

Package Size: #63631-1 50ul #63631-2 100ul


Description

Product Name	FGF2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FGF2 (AA: 189-288) expressed in E. Coli.
Target Name	FGF2
Other Names	BFGF; FGFB; FGF-2; HBGF-2
Accession No.	P09038
Calculated MW	30.8kd
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200-1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Hela cells using FGF2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.