

FGFR1 Mouse mAb

Catalog No: #63632

Package Size: #63632-1 50ul #63632-2 100ul


Description

Product Name	FGFR1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant extracellular fragment of human FGFR1 (aa22-376) fused with hlgGfc tag expressed in HEK293 cells.
Target Name	FGFR1
Other Names	CEK; FLG; OGD; FLT2; KAL2; BFGFR; CD331; FGFBR; HBGFR
Accession No.	P11362
Calculated MW	92kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ELISA:1/10000

Images



Western blot analysis using FGFR1 mouse mAb against extracellular domain of human FGFR1 (aa22-376).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.