

FGFR3 Mouse mAb

Catalog No: #63633

Package Size: #63633-1 50ul #63633-2 100ul


Description

Product Name	FGFR3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FGFR3 (AA: 529-694) expressed in E. Coli.
Target Name	FGFR3
Other Names	ACH; CEK2; JTK4; CD333; HSFGR3EX
Accession No.	P22607
Calculated MW	87.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A549 cells using FGFR3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using FGFR3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using FGFR3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.