

FGG Mouse mAb

Catalog No: #63636

Package Size: #63636-1 50ul #63636-2 100ul


Description

Product Name	FGG Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FGG expressed in E. Coli.
Target Name	FGG
Other Names	FGG
Accession No.	P02679
Calculated MW	52kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 FC:1/200 - 1/400 ELISA:1/10000


Images




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Immunofluorescence analysis of 3T3-L1 cells using FGG mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of HepG2 cells using FGg mouse mAb (blue) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using FGg mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using FGg mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.